



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® compact shaker KS 15 B control

€2,926.00
plus VAT &
Shipping

Product Images



Description

Description For reproducible work Scope of delivery Basic unit without vibrating plate, without attachment With touch display and programmable control Technical data Shaking motion: back and forth (horizontal) Shaking surface: 400 x 300 mm Max. Load: 15 kg Speed control range: 30 - 420 min⁻¹, in steps of 5 Stroke: 17 mm Running time: programmable / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: approx. 10 - 15 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 510 x 490 x 150 mm Weight: 19 kg Loading possibilities (examples) Erlenmeyer flasks 100 ml 20 pieces Erlenmeyer flasks 250 ml 15 pieces Erlenmeyer flasks 500 ml 12 pieces Erlenmeyer flasks 1000 ml 5 pieces Separating funnel 250 ml 1 piece

Additional Information

No.	NB-0142
Manufacturer (Brand)	Bühler
VGKL number	526161720
Transport temperature	Room temperature

