



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® Universal Shaker SM 30 B control

€4,278.00
plus VAT &
Shipping

Product Images



Description

Description Benchtop shaker with freely programmable sequence for mixing difficult media Scope of delivery Basic unit including vibrating plate and rubber mat, without attachment frame With touch display and programmable control system Technical data Shaking motion: back and forth (horizontal) Shaking surface: 560 x 400 mm Max. Load: 30 kg Speed control range: 15 - 300 min⁻¹, in steps of 5 Stroke: 30 mm 46 mm (option: No. 0052 212) 50 mm (option: No. 0052 216) Running time: programmable / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: approx. 20 - 30 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % Outer dimensions (W x D x H): 680 x 610 x 160 mm Weight: 33 kg Loading possibilities (examples) Erlenmeyer flasks 100 ml 31 pieces Erlenmeyer flasks 250 / 300 ml 24 pieces Separating funnels 250 ml 6 pieces Separating funnels 1000 ml 2 pieces Test tube racks 5 pieces

Additional Information

No.	NB-0118
Manufacturer (Brand)	Bühler
VGKL number	526161020
Transport temperature	Room temperature

