



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® multi-piston shaker VKS 75 B control

€7,666.00
plus VAT &
Shipping

Product Images



Description

Description Robust large shaker for high loads in long-term operation (to be set up on the floor only) Scope of delivery Basic unit without vibrating plate, without attachment frame With touch display and programmable control system Technical data Shaking motion: back and forth (horizontal) Shaking surface: 760 x 600 mm Max. Load: 75 kg Speed control range: 20 - 200 min⁻¹, in steps of 5 Stroke: 50 mm 36 mm (option: No. 0052 213) 60 mm (option: No. 0052 214) 80 mm (option: No. 0052 215) Running time: programmable / continuous operation Electrical connection: 100 - 240 V, 50 - 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: max. 45 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 1050 x 835 x 250 mm Weight: 110 kg Loading possibilities (examples) Erlenmeyer flasks 500 ml 34 pieces Erlenmeyer flasks 2000 ml 12 pieces Separating funnels 1000 ml 6 pieces

Additional Information

No.	NB-0125
Manufacturer (Brand)	Bühler
VGKL number	526161120
Transport temperature	Room temperature

