



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lauda® Varioshake VS 8 BE Shaker 230 V, 50/60 Hz

€1,240.00
plus VAT &
Shipping

Product Images



Description

LAUDA Varioshake VS 8 BE Shaker 230 V; 50/60 Hz Performance features: Shaker with bidirectional movement, high load capacity, low noise in compact design Electronic speed control, infinitely adjustable, smooth start-up Clear instrument panel for easy operation with rotary knobs Timer adjustable up to 60 min. or continuous operation Movement table made of anodized aluminum with four plastic pins Drive by an overload-protected AC motor Compact, low-wear mechanics with mass balancing for safe standing Constant speed, load-independent Excellent for continuous operation Outer housing made of electrolytically galvanized sheet steel, powder-coated Technical data: Ambient temperature range: 10 ... 60 °C Motion table (WxD): 330 x 330 mm Power consumption max.: 0.07 kW Shaking amplitude: 20 mm Shaking frequency: 20 ... 300 1/min Load capacity max.: 8 kg Dimensions (WxDxH): 350 x 375 x 160 mm Weight: 11 kg Power plug: Power cord with angled Schuko plug (CEE7/7) Power supply: 230 V; 50/60 Hz

Additional Information

No.	LA-1277
Manufacturer (Brand)	Lauda
VGKL number	101103056
Transport temperature	Room temperature

