



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® Microplate Shaker TiMix 2

€1,793.00
plus VAT &
Shipping

Product Images



Description

Description Especially for small sample volumes in MT plates Scope of delivery Basic unit Basic unit including attachment for 4 microtiter plates Technical data Shaking motion: exact circular Shaking surface: 200 x 295 mm for 4 MT plates Max. Loading: 2 kg Speed control range: 100 - 1400 min⁻¹, infinitely variable Stroke: 3 mm Timer: 0 - 120 min / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz please specify when ordering Protection class: IP 21 Heat output: approx. 15 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 220 x 410 x 110 mm Weight: 8.5 kg Loading options (examples) 4 microtiter plates Deepwell or other plates in MT plate format

Additional Information

No.	NB-0123
Manufacturer (Brand)	Bühler
VGKL number	526161100
Transport temperature	Room temperature

