



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Büehler® microplate shaker TiMix 5 control

€3,464.00
plus VAT &
Shipping

Product Images



Description

Description Top technology - also for use in robotic analysis systems Scope of delivery Basic unit without oscillating plate, without attachment With touch display and programmable control system Technical data Shaking motion: exact circular Shaking surface: 400 x 300 mm for 8 MT plates Max. Load: 5 kg Speed control range: 100 - 1400 min⁻¹, in steps of 5 Stroke: 3 mm (special stroke on request) Running time: programmable / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Power consumption: max. 30 W Heat output: approx. 11 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 510 x 490 x 150 mm Weight: 20 kg Loading options (examples) Standard attachment 0052 101 8 MT plates Attachment with pins 0052 096 8 MT plates Standard attachment 0052 101 + extension tray 0052 102 16 MT plates Attachment with maxi pins 0052 104 24 MT plates Attachment with metal pins 0052 125 48 MT plates

Additional Information

No.	NB-0139
Manufacturer (Brand)	Bühler
VGKL number	526161670
Transport temperature	Room temperature

