




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bühler® Universal Shaker SM 30 C

€4,109.00
plus VAT &
Shipping

Product Images



Description

Description Bench shaker for high loads in continuous operation Technical data Shaking motion: switchable circular or back + forth Shaking surface: 560 x 400 mm Max. Load: 30 kg Speed control range: 15 - 300 min⁻¹, infinitely variable Stroke: 26 mm Timer: 0 - 120 min / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: approx. 20 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 680 x 610 x 160 mm Weight: 33 kg Loading possibilities (examples) Erlenmeyer flasks 100 ml 42 pieces Erlenmeyer flasks 250 / 300 ml 20 pieces Separating funnels 250 ml 6 pieces Separating funnels 1000 ml 4 pieces Test tube racks 5 pieces

Additional Information

No.	NB-0122
Manufacturer (Brand)	Bühler
VGKL number	526161050
Transport temperature	Room temperature

