




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bühler® Universal Shaker SM 30 C control

**€5,222.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Description Benchtop shaker with freely programmable sequence for mixing difficult media Scope of delivery Basic unit including vibrating plate and rubber mat, without attachment frame With touch display and programmable control system Technical data Shaking motion: switchable circular or back + forth Shaking surface: 560 x 400 mm Max. Load: 30 kg Speed control range: 15 - 300 min<sup>-1</sup>, in steps of 5 Stroke: 26 mm Running time: programmable / continuous operation Electrical connection: 230 V or 115 V, 50 / 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: approx. 20 - 30 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 680 x 610 x 160 mm Weight: 33 kg Loading possibilities (examples) Erlenmeyer flasks 100 ml 31 pieces Erlenmeyer flasks 250 / 300 ml 24 pieces Separating funnels 250 ml 6 pieces Separating funnels 1000 ml 2 pieces Test tube racks 5 pieces

## Additional Information

---

No.	NB-0121
Manufacturer (Brand)	Bühler
VGKL number	526161040
Transport temperature	Room temperature

