



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bühler® multi-piston shaker VKS 75 A control

**€7,538.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Description Robust large shaker for high loads in long-term operation (to be set up on the floor only) Scope of delivery Basic unit without vibrating plate, without attachment frame With touch display and programmable control system  
Technical data Shaking motion: circular Shaking surface: 760 x 600 mm Max. Load: 75 kg Speed control range: 20 - 200 min<sup>-1</sup>, in steps of 5 Stroke: 26 mm Running time: programmable / continuous operation Electrical connection: 100 - 240 V, 50 - 60 Hz (please specify when ordering) Protection class: IP 21 Heat output: max. 45 W Ambient temperature: 5°C to 50°C Permissible humidity: ~ 85 % External dimensions (W x D x H): 1050 x 835 x 250 mm Weight: 110 kg Loading possibilities (examples) Erlenmeyer flasks 500 ml 34 pieces Erlenmeyer flasks 2000 ml 12 pieces Separating funnels 1000 ml 6 pieces

## Additional Information

---

No.	NB-0124
Manufacturer (Brand)	Bühler
VGKL number	526161110
Transport temperature	Room temperature

