



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Accessories for MaxQ 7000 Water Bath Shaker Movable dolly for underLab

€1,040.00
plus VAT &
Shipping

Product Images



Description

Accessories for MaxQ 7000 water bath shakers Thermo Scientific MaxQ 7000 water bath shaker accessories offer optimal flexibility. With exceptional innovation, reliability and accuracy, Thermo Scientific MaxQ 7000 shakers, platforms and accessories are the optimal choice for your laboratory. Thermo Scientific MaxQ 7000 water bath shaker accessories offer optimal flexibility. With exceptional innovation, reliability and accuracy, Thermo Scientific MaxQ 7000 shakers, platforms and accessories are the optimal choice for your laboratory. Description: Movable dolly for underLab bench operation Item description: Movable dolly for underLab bench operation For use with: For MaxQ 7000 shakers Type: Movable Dolly For use with (devices): MaxQ 7000

Additional Information

No.	LD-5025
Manufacturer (Brand)	Thermo Scientific

