



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific Digital Shake Dry Bath

**€1,530.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Digital Shaker Dry Bath A Thermo Scientific Digital Heated Shaker Dry Bath is suitable for use as a dry bath, shaker, or both. Connector type: cable set with various connectors Certification/Compliance: CE Controller type: control with PID Diameter (metric) Orbit: 2 mm Capacity (metric): 0.5 ml; 1.5 ml, 50 ml Dimensions (L x W x H): 13.1 x 10.9 x 6.7 in. (334 x 276 x 170 mm) Temperature range (metric): Ambient temperature: +5 °C to 100 °C Display: digital Electrical requirements: 230 V, 50/60 Hz Frequency: 50/60 Hz Number of blocks: Not supplied, purchase separately Shipping Weight (metric): 8.3 kg Shipping Weight (Anglo): 18.3 lb. Speed: 150 to 1,500 rpm (96-well block); 150 to 1,200 rpm (blocks for 0.5 ml, 1.5 ml, and 50 ml tubes) Temperature Accuracy: 0.5 °C (0.9 °F) at 80 °C Warranty: 2 years Illumination: LED Type: Platform Shaker No. of Places: 96 Unit Size: Each

## Additional Information

---

No.	LD-2716
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

