



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific Teleshake Magnetic Shaking System

**€1,980.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Teleshake Magnetic Shaking System Take advantage of powerful shaking applications at 100 to 2000 rpm with this remote-controlled, wear-free magnetic shaking system. Height (Anglo): 1.85 in. Height (metric): 47 mm Width (anglo): 5.74 in. Width (metric): 146 mm Certification/Conformity: certificate of analysis (COA), certificate of conformity (COC) Type: shaker Amplitude: 2 mm Capacity: 1 Item description: control unit (230 V EU) included Humidity during operation: +5 to +40 °C at 95% RH Speed: 100 to 2,000 rpm Voltage: 230 V EU Temperature range (metric): 5 deg;C to 40 deg;C (operating temperature) Unit Size: Each

## Additional Information

---

No.	LD-1147
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

