



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific stainless steel ion exchange cartridges (pressure resistant)

€1,170.00
plus VAT &
Shipping

Product Images

Description

Stainless steel ion exchange cartridges (pressure resistant) These stainless steel cartridges are proven mixed bed systems for demineralized water with flow rates from 300 to 3000 l/h and are supplied ready for connection.
Description: DI 2800 Ion Exchange Stainless Steel (pressure resistant) Cartridges Purity: 0.1 to 20 µS/cm Type: Ion Exchange Cartridge Capacity: 2,800 L Flow Rate: 950 l/h Purification Method: Ion Exchange Resin Type: Mixed Bed Unit Size: Each

Additional Information

No.	LD-0035
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

