




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Stakpure® pure and ultrapure water system OmniaLab ED+ 20

€12,640.00
plus VAT &
Shipping

Product Images



Description

OmniaLab ED+ 20 pure and ultrapure water system for the complete laboratory water supply of autoclaves, dishwashers and ultrapure water systems by continuous self-regenerating electrodeionization, mobile as tower unit or space-saving in a laboratory base cabinet, standard incl. 100 l pure water tank with recirculation and pressure outlet, OptiFill dispenser & monitoring unit and leakage sensor Pure water ASTM type II: Conductivity $\mu\text{S}/\text{cm}$: 0.067-1 Resistance $\text{M}\Omega \times \text{cm}$: 15-1 Pure water capacity l/h: 20 Ultrapure water ASTM type I: Conductivity $\mu\text{S}/\text{cm}$: 0.055 Resistance $\text{M}\Omega \times \text{cm}$: 18.2 TOC value ppb: 5-10 Bacteria content KbE/ml: < 0.1 Dimensions with dispenser mm: W 511 x D 575 x H 1520 Dimensions undersink unit tank mm: W 511 x D 575 x H 800

Additional Information

No.	PE-0024
Manufacturer (Brand)	Stakpure
VGKL number	104541872
Transport temperature	Room temperature

