




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Stakpure® Pure and Ultrapure Water System OmniaTap 12 UV/UF

€7,840.00
plus VAT &
Shipping

Product Images



Description

OmniaTap 12 UV/UF pure and ultrapure water system for direct connection to the drinking water line for the production of pure and ultrapure water according to ASTM Type II and I as standard incl. 10 l pure water tank with recirculation and optional pressure outlet, OptiFill dispenser & monitoring unit and leakage sensor, completely ready to connect incl. Pretreatment cartridge, ultrapure water cartridge, UV unit 185/254 nm, ultrafiltration and sterile filter capsule 0.2 µm ultrapure water ASTM type I: Conductivity µS/cm: 0.055 Resistance MΩ x cm: 18.2 TOC value ppb: 1-5 Endotoxins (pyrogens) EU/ml: 0.001 Bacteria content KbE/ml: < 0.1 Withdrawal rate l/min: up to 2 Pure water ASTM type II: Conductivity µS/cm: 0.067-0.1 Resistance MΩ x cm: 15-10 Pure water rate l/h: 12 Dimensions with dispenser mm: W 390 x D 615 x H 720

Additional Information

No.	PE-0023
Manufacturer (Brand)	Stakpure
VGKL number	104548213
Transport temperature	Room temperature

