



Thermo Scientific Barnstead™ GenPure™

Product Images

Description

Barnstead™ GenPure™ The Thermo Scientific™ Barnstead™ GenPure™ system produces type 1 ultrapure water of consistent quality for your most demanding and sensitive applications. Frequency: 50/60 Hz Feedwater Connector: R 3/4 in. Type of Product Water: Type 1 TOC Monitor with UV Intensity Monitoring: No Type: Water System Electrical Resistance: 18.2 Mohm.cm Conductivity: 0.055 S / cm Feedwater Source: Prepared Flow Rate: Up to 2 l/min TOC (Total Organic Carbon): 1 to 5 ppb Endotoxin Concentration: 0.001 EU/mL Bacterial Content: It; 1CFU/mL Particles 0.2 µm/mL: It;1 Feedwater Monitoring: yes Feedwater Pressure: 1 to 6 bar Operating Pressure: 2 to 6 bar Voltage: 90/240 V Temperature Range (metric): 2 °C° to 35 °C° Description: GenPure ultrapure water system with UV-photo-oxidation and ultrafiltration membrane Power consumption: 0.1 kW Electrical requirements: 90/240 V 50/60 Hz Unit Size: Each

Additional Information

No.	LD-1253
Manufacturer (Brand)	Thermo Scientific

