



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific Barnstead™ GenPure™ xCAD Plus

**€11,800.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Barnstead™ GenPure™ xCAD Plus Dispense Type 1 ultrapure water from up to three dispensers simultaneously with the Thermo Scientific Barnstead GenPure xCAD Plus high purity water treatment system. Frequency: 50/60 Hz Feedwater Connector: 3/4 in. Type of Product Water: Type 1 TOC Monitor with UV Intensity Monitoring: Yes Type: Water Treatment System Electrical Resistance: 18.2 Mohm.cm at 25 °C Conductivity: 0.055µS/cm Flow Rate: Up to 2 L/min TOC (Total Organic Carbon): 1 to 5 ppb Bacterial Content: It;0.01 CFU/mL Particles 0.2 µm/ml: It;1 Remote Dispenser: Yes Feedwater Monitoring: Integrated Height (Anglo-American): 24 in. Width (anglo): 14.6 in. Operating Pressure: 0.1 to 6 bar Voltage: 90/240 V Temperature Range (metric): 2 °C° to 35 °C° Depth (Anglo-American): 13 in. Description: Barnstead GenPure xCAD Plus UV - TOC (wall version) Power Consumption: 0.1 kW Electrical Requirements: 90/240 V 50/60 Hz Unit Size: Each

## Additional Information

---

No.	LD-1301
Manufacturer (Brand)	Thermo Scientific
VGKL number	108662226
Transport temperature	Room temperature

