



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Barnstead™ LabTower™ EDI

€18,200.00
plus VAT &
Shipping

Product Images

Description

Barnstead™ LabTower™ EDI Safe and convenient storage of pure water with the Barnstead LabTower EDI water treatment system with integrated 100 l tank. Frequency: 50/60 Hz Feedwater Connector: 3/4 in. NPT Type of Product Water: Type 1 and 2 TOC Monitor with UV Intensity Monitoring: No Type: Water Treatment System Electrical Resistivity: Type 1: 18.2 Mohm.cm at 25 °C; Type 2: 10 to 15 Mohm.cm at 25 °C Conductivity: Type 1: 0.055 µS/cm; Type 2: 0.067 to 0.1 µS/cm Feed water source: Pipe Flow rate: 1.5 l/min TOC (Total Organic Carbon): 1 to 5 ppb Bacterial content: It; 0.01 CFU/ml Particles 0.2 µm/ml: It; 1 Feed water monitoring: Yes Feed water pressure: 1 to 6 bar Operating pressure: 2 to 6 bar Voltage: 90/240 V Bacterial retention: 99 % Description: LabTower EDI 15 l/h with integrated 100 l tank Permeate capacity: 15 l/h at 15 °C Power consumption: 0.25 kW Electrical requirements: 90/240 V 50/60 Hz Unit Size: Each

Additional Information

No.	LD-1274
Manufacturer (Brand)	Thermo Scientific
VGKL number	108664415
Transport temperature	Room temperature

