



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Barnstead™ Pacific™ RO

€6,510.00
plus VAT &
Shipping

Product Images

Description

Barnstead™ Pacific™ RO Functional and economical production of high purity water with a modular design. Bacterial retention: 99% Description: 0.5 ; 1.0; 1.5 ; 2.0; 4.0nm Electrical requirements: 230 V, 50/60 Hz Feed water connection: 3/4 in. NPT Frequency: 50/60 Hz Operating pressure: 1.5 to 6 bar Permeate capacity: 12 l/h Power consumption: 0.1 kW Temperature range (metric): 2 °C° to 35 °C° Feed water monitoring: Integrated Flow rate: 12 l/h Type: Reverse osmosis clean water system Unit Size: Each

Additional Information

No.	LD-1271
Manufacturer (Brand)	Thermo Scientific
VGKL number	108662012
Transport temperature	Room temperature

