



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Stakpure® Pure and Ultrapure Water System OmniaTap-T xs touch 8 UV-TOC

€7,720.00
plus VAT &
Shipping

Product Images

Description

OmniaTap-T xs touch 8 UV-TOC pure and ultrapure water system for direct connection to the drinking water line for the production of pure and ultrapure water according to ASTM Type II and I as standard incl. 7 l pure water tank with recirculation and optional pressure outlet, OptiFilltouch dispenser & monitoring unit and leakage sensor, completely ready to connect incl. Pretreatment cartridge, ultrapure water cartridge, UV unit 185/254 nm and sterile filter capsule 0.2 µm ultrapure water ASTM type I: Conductivity µS/cm: 0.055 Resistance MΩ x cm: 18.2 TOC value ppb:

Additional Information

No.	PE-0136
Manufacturer (Brand)	Stakpure
Transport temperature	Room temperature

