



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Thermo Scientific Barnstead™ LabTower™ RO

**€14,200.00**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

---

Barnstead™ LabTower™ RO Convert tap water to reverse osmosis (RO) water with the Thermo Scientific™ Barnstead™ LabTower™ RO system with integrated 100L storage tank. Bacterial content: 99% Description: LabTower RO 40L/hr. Electrical requirements: 100/240 V 50/60 Hz Feed water connection: 3/4 in. NPT Feed water pressure: 1 to 6 bar Frequency: 50/60 Hz Operating pressure: 2 to 6 bar Particle 0.2 µm/mL: 1 µg/mL Permeate capacity: 40 l/h at 10 °C Power consumption: 0.25 kW Silica removal: 98% Temperature range (metric): 2 °C to 35 °C Feed water monitoring: Yes Feed water source: Cone Flow rate: 40 l/h at 10 °C TOC Monitor with UV Intensity Monitoring: No Type: Reverse osmosis clean water system Unit Size: Each

## Additional Information

---

No.	LD-1273
Manufacturer (Brand)	Thermo Scientific
VGKL number	108664440
Transport temperature	Room temperature

