



Thermo Scientific Barnstead™ Pacific™ RO

Product Images

Description

Barnstead™ Pacific™ RO Functional and economical production of high purity water with a modular design. Bacterial retention: 99%
Description: X3FR (cooling); 230V Electrical requirements: 230 V, 50/60 Hz Feed water connection: 3/4 in. NPT Frequency: 50/60 Hz
Operating pressure: 1.5 to 6 bar Permeate capacity: 7 l/h Power consumption: 0.1 kW Temperature range (metric): 2 °C° to 35 °C°
Feed water monitoring: Integrated Flow rate: 7 l/h Type: Reverse osmosis clean water system Unit Size: Each

Additional Information

No.	LD-1270
Manufacturer (Brand)	Thermo Scientific
VGKL number	108662007

