



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Memmert® Water Bath WNB7, 7l, 10-95°C

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Water bath WNB7, 7 l, working temperature range min. 5 °C above ambient temperature up to +95 °C and additional boiling step, setting temperature range +10 °C to +95 °C; electronic microprocessor PID controller with self-diagnostic system, 1 Pt100 sensor, integrated digital clock from 1 min to 99.59 h for: On (continuous operation) / Wait (delayed On for continuous and timed operation) / Hold, 2-fold overtemperature protection, usable space and housing made of stainless steel; usable space: 240 x 210 x 140 mm, housing: 468 x 356 x 238 mm, voltage: 230 V, 50/60 Hz, max. power consumption: approx. 1200 W

Additional Information

| | |
|-----------------------|------------------|
| No. | MM-11054 |
| Manufacturer (Brand) | Memmert |
| Transport temperature | Room temperature |

