



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Clamps and Fasteners for Precision™ Water Baths Each NA Clip Fastener

€124.00
plus VAT &
Shipping



Description

Clamps and Attachments for Precision™ Water Baths Clamps and attachments for use with Thermo Scientific™ Precision™ water baths.™ Bottle Clamps Secure flasks of various sizes to bath platform Stainless Steel Each requires a fixture (Catalog # 3166189 - order separately) Tube Clamp Stainless Steel For 13 to 25 mm tubes Each requires a fixture (Catalog # 3166189 - order separately) Custom Group: Miscellaneous Water Bath Accessories CustomPDPBanner: LegalEntityNotice LeadTargetGroup: xx_ELMS Product Size: Each Search Display: Family Umbrella Brand: Thermo Scientifictrade; Capacity: NA For use with: For use in Precision Reciprocal, Dubnoff, and Shallow-Form Shaking Water Baths. Type: clip closure For use with (devices): Precision Reciprocal Dubnofftrade; Shaking Water Bath For use with (devices): Precision Reciprocal Shallow-Form Shaking Water Bath Custom Group: Miscellaneous Water Bath Accessories CustomPDPBanner: LegalEntityNotice LeadTargetGroup: xx_ELMS Product Size: Each Search Display: Family Umbrella Brand: Thermo Scientifictrade; Capacity: NA For Use With: For use in Precision Reciprocal, Dubnoff, and Shallow-Form Shaking Water Baths. Type: Clip Fastener For Use With (Equipment): Precision Reciprocal Dubnofftrade; Shaking Water Bath For Use With (Equipment): Precision Reciprocal Shallow-Form Shaking Water Bath

Additional Information

No.	LD-0513
Manufacturer (Brand)	Thermo Scientific
Transport temperature	Room temperature

