



Thermo Scientific AC 200-S7 230V

Product Images

Description

AC 200-S7 (230V/50Hz) HEATING BATH CIRCULATOR STEEL, AC200-S7, Amb+13 TO +200 C, 4-7 L BATH VOLUME, 2.0 kW HEATER, FORCE SUCTION PUMP 20 L/MIN REMOTE SENSOR PORT; USB; RAMP PROGRAM; MULTIFUNCTION PORT; 8 LANGUAGES; UL/CSA,

Additional Information

No.	LD-4276
Manufacturer (Brand)	Thermo Scientific
VGKL number	105336107

