



Thermo Scientific SC150-A10B 230V

Product Images

Description

SC 150-A10B (230V/50Hz) OPEN REFRIGERATION BATH CIRCULATOR STEEL, SC150-A10B, -10 to +100 DEGREE C, 250 W COOLING CAPACITY AT 20C, 17-30 L BATH VOLUME, 2.0 kW HEATER, PUMP 17 L/MIN USB;

Additional Information

No.	LD-4244
Manufacturer (Brand)	Thermo Scientific
VGKL number	105333410

