



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Krüss® Density Meter – Set 2 with glass bending oscillator, with syringe in the

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

Density meter - Set 2 with glass bending oscillator, with syringe into density meter DS7700 - including drying unit with 2/2-way valve and small Tygon and PTFE tube sets Measuring range: 0-3 g/cm<sup>3</sup> Accuracy: +/-0.001 g/cm<sup>3</sup> Temperature range: 10-40 °C Temperature stability: +/-0.02 °C Sample volume: approx. 0.9 ml Air pressure sensor: 800-1100 hPa Accuracy: +/-2 hPa Interfaces: RS-232, USB & Ethernet Power supply: 100-240 V, 50/60 Hz, 60 W Dimensions (WxHxD): 22x22x43 cm<sup>3</sup>

## Additional Information

---

No.	KS-0515
Manufacturer (Brand)	Krüß
VGKL number	106117702
Transport temperature	Room temperature

