



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Krüss® Density Meter - Set 3 with Glass Bending Transducer, Semi-Automatic Sample Feeder

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Density meter - Set 3 with glass bending transducer, Semi-automatic sample feeding with peristaltic pump into the density meter DS7800 - including drying unit with 3/2-way valve and large Tygon tubing set Measuring range: 0-3 g/cm³ Accuracy: +/-0.0001 g/cm³ Temperature range: 10-40 °C Temperature stability: +/-0.02 °C Sample volume: approx. 0.9 ml Air pressure sensor: 800-1100 hPa Accuracy: +/-2 hPa Interfaces: RS-232, USB & Ethernet Power supply: 100-240 V, 50/60 Hz, 60 W Dimensions (WxHxD): 22x22x43 cm³

Additional Information

No.	KS-0521
Manufacturer (Brand)	Krüss
VGKL number	106117803
Transport temperature	Room temperature

