



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Krüss® Density Meter – Set 4 with Glass Bending Oscillator, Fully Automatic Sample Supply

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

Density meter - Set 4 with glass bending transducer, Fully automatic sample feeding with peristaltic pump into the density meter DS7700 - including drying unit with 3/2-way valve, autosampler and large Tygon tubing set Measuring range: 0-3 g/cm<sup>3</sup> Accuracy: +/-0.001 g/cm<sup>3</sup> Temperature range: 10-40 °C Temperature stability: +/-0.02 °C Sample volume: approx. 0.9 ml Air pressure sensor: 800-1100 hPa Accuracy: +/-2 hPa Interfaces: RS-232, USB & Ethernet Power supply: 100-240 V, 50/60 Hz, 60 W Dimensions (WxHxD): 22x22x43 cm<sup>3</sup>

## Additional Information

---

No.	KS-0517
Manufacturer (Brand)	Krüß
VGKL number	106117704
Transport temperature	Room temperature

