



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Krüss® Density Meter - Set 5 with Glass Bending Oscillator, Fully Automatic Sample Loading

Price on
request

plus VAT &
Shipping

Product Images



Description

Density meter - Set 5 with glass bending oscillator, Fully automatic sample feeding with peristaltic pump into the density meter DS7700 - including drying unit with 3/2-way valve, autosampler and large PTFE tubing set Measuring range: 0-3 g/cm³ Accuracy: +/-0.001 g/cm³ Temperature range: 10-40 °C Temperature stability: +/-0.02 °C Sample volume: approx. 0.9 ml Air pressure sensor: 800-1100 hPa Accuracy: +/-2 hPa Interfaces: RS-232, USB & Ethernet Power supply: 100-240 V, 50/60 Hz, 60 W Dimensions (WxHxD): 22x22x43 cm³

Additional Information

No.	KS-0518
Manufacturer (Brand)	Krüß
VGKL number	106117705
Transport temperature	Room temperature

