



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Sartorius® Infrared Moisture Meter MA35

€2,880.00
plus VAT &
Shipping

Product Images

Description

The MA35 is the ideal entry-level model for fast and reliable determination of the material moisture of liquid, pasty and solid substances using the thermogravimetry method. The selection of durable, robust components, insensitive to dirt and vibration, guarantees reliable measurement results even in harsh environments. An ideal partner for use in production, incoming goods or in the laboratory. Advantages Time-saving alternative to the drying oven - measurement results in just 10 - 20 minutes Simple operation Robust design Ideal for simple routine checks End-point determination of moisture measurement using fully automatic or time switch-off Best price/performance ratio in Sartorius quality - Made in Germany

Additional Information

No.	ST-15538
Manufacturer (Brand)	Sartorius
VGKL number	105551031
Transport temperature	Room temperature

