



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Sartorius® Infrared Moisture Meter, 70g, 1mg

**€4,320.00**  
**plus VAT &**  
**Shipping**

### Product Images

## Description

---

The MA37 is the new Sartorius basic model for moisture analysis. It works faster and is easier to operate than its predecessor, the MA35. In addition, it has a compact design and high measurement accuracy. Its BetterClean design also allows easy cleaning and the status light indicates the current status of the measurement. The MA37 is therefore ideal for operations characterized by consistent measurement tasks and guides the user step-by-step through the measurement process. The MA37 moisture analyzer is used for fast and reliable determination of the material moisture of liquid, pasty and solid substances using the thermogravimetric method. For daily recurring moisture determination tasks such as in production or incoming goods inspection, the MA37 is perfectly suited. Advantages: High measurement speed Effortless cleaning User-friendly operation Display of measurement progress Downward compatibility with the previous model MA35

## Additional Information

---

No.	ST-15539
Manufacturer (Brand)	Sartorius
Transport temperature	Room temperature

