



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Testo® 608-H2 - Thermohygrometer

**€126.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Gardening; offices; warehouses, laboratories, living rooms: Wherever permanent but discreet room climate monitoring is required, the testo 608-H2 thermohygrometer is ideal for use. Thanks to the hanging and positioning device, it can be mounted on the wall or placed on tables and shelves. testo 608-H2: the practical thermohygrometer from Testo

The testo 608-H2 thermohygrometer is the right choice if you do not need extensive documentation and archiving of measurement data, but simply want to permanently measure and view the current room climate values: Temperature, humidity and dew point are shown on the large display, and the values are easy to read even at a distance. In addition to the max. and min. value display, the testo 608-H2 has an optical LED alarm that reports when limit values are exceeded. Both upper and lower limit values can be set individually. The precise and long-term stable humidity sensor makes the testo 608-H2 a reliable; low-maintenance measuring instrument; which can still be relied on after years. The thermohygrometer has a battery life of approximately one year. For safety, the testo 608-H2 has integrated battery monitoring; a critical battery status is shown on the display in good time.

## Additional Information

---

No.	TO-0096
Manufacturer (Brand)	Testo
VGKL number	100946082
Transport temperature	Room temperature

