




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell® replacement pocket thermometer, enclosed type, -35+50:1°C, capillary prismatic

€25.79 As low
as
€11.25
plus VAT &
Shipping



Description

Spare pocket thermometer, enclosed type, -35+50:1°C, capillary prismatic uncoated, red special filling, 140x8,5mm, fully immersed

Additional Information

No.	LL-2241
Manufacturer (Brand)	Amarell

