




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell® cold laboratory thermometer, DIN 12778, enclosed type, -100+30:1°C, capillary

€57.31 As low as
€41.67
plus VAT &
Shipping



Description

Cold laboratory thermometer, DIN 12778, enclosed type, $-100+30:1^{\circ}\text{C}$, capillary prismatic uncoated, toluene filling, 300x8-9mm, fully immersing, Richter lock with knob, suitable for official certification

Additional Information

No.	LL-2538
Manufacturer (Brand)	Amarell
VGKL number	161019006

