



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell® cold laboratory thermometer, DIN 12778, rod type, -50+50:1°C, white coated


€44.30 As low

as

€32.32

**plus VAT &
Shipping**

Product Images



Description

Cold laboratory thermometer, DIN 12778, stem type, -50+50:1°C, white coated, toluene filling, 300x6-7mm, full immersion, top with ring, suitable for official certification

Additional Information

| | |
|----------------------|---------|
| No. | LL-2555 |
| Manufacturer (Brand) | Amarell |

