



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell® laboratory thermometer with saddle, enclosed type, -38+50:0.1°C

**Price on
request
plus VAT &
Shipping**

Amarell

Precision-Thermometers and Hydrometers

- General Purpose Thermometers
- Precision Thermometers
- Precision Thermometers in Sets
- Industrial Thermometers
- Cylinder Thermometers
- Dial Thermometers
- Quartz Glass Thermometers
- Thermometers for Calorimetry
- Thermometers for Microbiology
- Thermometers for Special Purposes
- Thermometers with Standard Ground Joint
- Thermometers According to DIN, ASTM, IP, BS, SFTTC and AFNOR
- AMB-Therm[®] - non-Wetting Low-Range Precision Thermometers
- AMB-Galvan[®] - non-Wetting Low-Range Precision Thermometers

- General Purpose Hydrometers
- Hydrometers According to National and International Standards
- Hydrometers for Mineral Oil and Liquefied Gas Testing
- Hydrometers for BSA, Water According to Baumé and for Sugar Solutions (Saccharometers)
- Density Hydrometers and Alcoholometers
- Hydrometers for Special Applications
- Fittings for Hydrometric Measurements

PRECISION

Amarell GmbH & Co.KG Lindenstraße 43 D-87892 Kreuzwertheim Germany www.amarell.de

Description

Laboratory thermometer with saddle, enclosed type, $-38+50:0.1^{\circ}\text{C}$, with special capillary inner coating, with non-wetting blue special filling, 550x9-10mm, fully immersing, Richter closure with knob, with factory test certificate with 1 test point

Additional Information

No.	LL-3088
Manufacturer (Brand)	Amarell
Transport temperature	Room temperature

