



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## **Amarell® laboratory thermometer with saddle, enclosed type, 0+100:0.1°C**

**Price on  
request  
plus VAT &  
Shipping**

# Amarell

## Precision-Thermometers and Hydrometers

- General Purpose Thermometers
- Precision Thermometers
- Precision Thermometers in Sets
- Industrial Thermometers
- Cylinder Thermometers
- Dial Thermometers
- Quartz Glass Thermometers
- Thermometers for Calorimetry
- Thermometers for Microbiology
- Thermometers for Special Purposes
- Thermometers with Standard Ground Joint
- Thermometers According to DIN, ASTM, IP, BS, SFTTC and AFNOR
- AMB-Therm<sup>®</sup> - non-Wetting Low-Hazard Precision Thermometers
- AMB-Glass<sup>®</sup> - non-Wetting Low-Hazard Precision Thermometers

- General Purpose Hydrometers
- Hydrometers According to National and International Standards
- Hydrometers for Mineral Oil and Liquefied Gas Testing
- Hydrometers for BSA, Water According to Baumé and for Sugar Solutions (Saccharometers)
- Brewery Hydrometers and Alcohometers
- Hydrometers for Special Applications
- Filtrage for Hydrometric Measurements

PRECISION

Amarell GmbH & Co KG Lindenstraße 43 D-87892 Kreuzwertheim Germany [www.amarell.de](http://www.amarell.de)

## Description

---

Laboratory thermometer with saddle, enclosed type, 0+100:0.1°C, with special capillary inner coating, with non-wetting blue special filling, 580x9-10mm, fully immersing, Richter closure with button, with factory test certificate with 1 test point

## Additional Information

---

No.	LL-3091
Manufacturer (Brand)	Amarell
Transport temperature	Room temperature

