



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Amarell® Adjustable NS Thermometers, -10+150°C, capillary uncoated, red special filling

**Price on
request
plus VAT &
Shipping**

Amarell

Precision-Thermometers and Hydrometers

- General Purpose Thermometers
- Precision Thermometers
- Precision Thermometers in Sets
- Industrial Thermometers
- Cylinder Thermometers
- Dial Thermometers
- Quartz Glass Thermometers
- Thermometers for Calorimetry
- Thermometers for Microbiology
- Thermometers for Special Purposes
- Thermometers with Standard Ground Joint
- Thermometers According to DIN, ASTM, IP, BS, SFTTC and AFNOR
- AMB-Therm[®] - non-Wetting Low-Hazard Precision Thermometers
- AMB-Galium[®] - non-Wetting Low-Hazard Precision Thermometers

- General Purpose Hydrometers
- Hydrometers According to National and International Standards
- Hydrometers for Mineral Oil and Liquefied Gas Testing
- Hydrometers for BSA, Water According to Baumé and for Sugar Solutions (Saccharimeters)
- Density Hydrometers and Alcoholometers
- Hydrometers for Special Applications
- Fittings for Hydrometric Measurements

PRECISION

Amarell GmbH & Co KG Lindenstraße 43 D-87892 Kreuzwertheim Germany www.amarell.de

Description

Adjustable standard ground joint thermometer AMA-Exakt, enclosed type, -10+150:0,5°C, capillary prismatic uncoated, red special filling, upper part 300x10mm, lower part 300x8-8,5 mm, Richter closure without knob, with ground joint NS 14,5/23 and universal closure, calibrateable, with thread correction certificate

Additional Information

No.	LL-2948
Manufacturer (Brand)	Amarell
Transport temperature	Room temperature

