



Testo® 623 - Thermo hygrometer

Product Images



Description

Placed on the table or mounted on the wall: The thermo-hygrometer in elegant design measures temperature and humidity in living rooms; offices and commercial premises constantly and reliably. Thermo-hygrometer testo 623: technical features With the thermo-hygrometer testo 623, you are relying on a reliable monitoring monitor for the indoor climate. The intervals at which the device displays the measured values can be set by the user, this ranges from a minimum of 5 min to a maximum of 24 hours. The values of the last 90 days are automatically stored. With the integrated alarm function the reaching of limit values can be shown. The current values, date and time can be seen at any time on the large display, as well as the statistics of the last 90 days. Thus, the long-term analysis of the room climate is possible at any time without any problems - completely without data transmission and evaluation on the PC. The long-term stable sensor guarantees long-lasting quality - so that you can still rely on the reliability of your thermo-hygrometer in a few years. Whether the thermo-hygrometer stands on the table or the sideboard like a picture frame or whether it hangs on the wall; you decide: both is possible; mounting material is included in the scope of delivery.

Additional Information

No.	TO-0098
Manufacturer (Brand)	Testo
VGKL number	100940623

