



Testo® Mini Insertion Thermometer

Product Images



Description

Small insertion thermometer with great reliability: The mini thermometer from Testo is ideal for measurements in liquids - e.g. in the laboratory - or for air temperature measurements in heating and ventilation technology. Likewise, you can also use the mini thermometer as a food thermometer. With its insertion probe, it is excellent for measuring the core temperature of food products during food inspection. For all measurements, whether food control or laboratory sampling, the mini thermometer offers high accuracy and ease of use. Thus, the piercing thermometer becomes a practical all-rounder that you can always carry with you thanks to its small size. Practical features of the mini plunge thermometer On the digital display of the plunge thermometer, you can easily read the exact measured value (unlike analog food thermometers, whose pointer always bears the risk of reading errors). In addition, the display offers a battery level indicator and the possibility; to display the temperature in the unit °C or °F. The supplied protective tube for the insertion probe can also be used as a holder to attach the mini thermometer to your shirt pocket using the retaining clip. This allows you to conveniently take the insertion thermometer with you to any measuring location.

Additional Information

No.	TO-0070
Manufacturer (Brand)	Testo
VGKL number	100961110

