



Testo® Waterproof Mini Insertion Thermometer

Product Images



Description

Mini in price, maxi in performance: The waterproof mini insertion thermometer measures the core temperature of liquid, soft and powdery media just as reliably as the air temperature. After use, you can simply clean the insertion thermometer under running water, or even in the dishwasher if it is very dirty. It can then be stored safely in the protective tube supplied. The clip integrated into the protective tube holds the thermometer reliably in your shirt or jacket pocket. From checking roasts in the kitchen (do not leave in the oven!) to checking the temperature in production, nurseries, air conditioning units or the laboratory: Thanks to its ability to measure both the air temperature and the core temperature, the insertion thermometer is suitable for many different areas of application.

Additional Information

No.	TO-0072
Manufacturer (Brand)	Testo
VGKL number	100961113

