



Lovibond® SD305 pH/ORP (Set 2)

Product Images



Description

pH electrode type 226, temperature sensor Pt1000 Waterproof handheld meter for measuring pH/redox/temperature Portable handheld meter for determining pH values. It is reliable , waterproof (IP67) as well as robust, making it an ideal solution for on-site analysis. Technically advanced A particularly high resolution (0.001 pH / 0.1 mV) for maximum accuracy, automatic temperature compensation and buffer detection are further features of the SD305 pH/ORP. User-friendly operation Convenient operation is rounded off by a large, backlit display. This means that you have the measurement results directly in view in any situation. Data management Data collected by the data logger can be transferred to the computer via USB interface (serial or analog). Protective housing The protective armor guarantees safe operation at all times, even in extreme environmental conditions. The stand attached to the back of the housing can also be used as a hanger or belt clip. Additional equipment Status display for battery and electrode, battery change indicator, Auto-Off, Auto-Hold, Min./Max. value memory are also included as functions.

Additional Information

No.	CV-0269
Manufacturer (Brand)	Lovibond

