



Lovibond® SD315 Oxi (Set 3)

Product Images



Description

Oxygen sensor, Pt/Pb, galvanic, cable length: 30 m Waterproof handheld meter for measuring dissolved oxygen (O₂), O₂ concentration/O₂ saturation/temperature Portable handheld meter for determining dissolved oxygen. It is reliable, rugged and ideal for instant measurements without time-consuming polarization time thanks to the Clark oxygen sensor. No polarization The well-known Clark oxygen sensor enables instant measurement without polarization time. Versatile The SD315 Oxi is aimed at industry, laboratories and municipal facilities such as sewage treatment plants or waterworks. Data management A data logger and the corresponding logger software ensure simplified processing of the acquired data. To ensure that this works under any visibility conditions, the device is supplied with a backlit display. Protective housing The protective armor guarantees not only a secure grip but also protection against damage caused by falling. The stand attached to the back of the case can also be used as a hanger or belt clip. Alarm The integrated alarm warns you if the measured value should ever deviate from the predefined values. Auto-Hold The Auto-Hold function enables easy reading of the measured values.

Additional Information

No.	CV-0273
Manufacturer (Brand)	Lovibond

