



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



WLD-TEC® Inoculation loop holder, length 245mm, made of stainless steel

€61.80
plus VAT &
Shipping

Product Images



Description

Inoculation loop holder made of stainless steel, with union nut, fablich coded, for inoculation loop wire 0.6 - 1 mm diameter. Autoclavable Length: 245mm PU = 4 pieces

Additional Information

No.	WT-0065
Manufacturer (Brand)	WLD-TEC
Transport temperature	Room temperature

