



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Bürkle® Detectable scraper, blue, PS, sterile, 80 mm

**€33.15** As low as  
**€29.50**  
plus VAT &  
Shipping

### Product Images



## Description

The detectable SteriPlast® scraper made of blue polystyrene is produced in our class 7 (10,000) clean room, individually packaged for disposable use and then gamma sterilized. Stable and handy scraper for disposable use. The detectable scraper is suitable for scraping, smoothing or cleaning, as well as for removing residues. In addition to the blue color, these samplers have a special additive in the material that makes them visible in metal detectors or by X-rays. These products, although made of plastic, can therefore be sorted out by the usual systems for checking foreign bodies. PS detectable, blue Individually packaged and gamma-sterilized Production, assembly and packaging in accordance with cleanroom class 7 Complies with EU food and FDA guidelines Detectable by color scanner, metal detector or X-rays Can be used as part of HACCP/IFS/BRC foreign object management Specially for the food, feed and pharmaceutical industries

## Additional Information

No.	KL-0958
Manufacturer (Brand)	Bürkle
VGKL number	100458230
Transport temperature	Room temperature

