



## Bürkle® MiniSampler PTFE, complete

### Product Images



### Description

Ideal when the highest purity is required. MiniSampler PTFE complete with accessories in carrying case, contents: vacuum pump, 1 sample bottle PFA 180 ml, 5 m PTFE suction hose 6 x 8 mm, 1 hose cutter, 1 stainless steel hose weight made of V2A (1.4301)/nickel-plated brass, 20 seal-it label seals. With the vacuum sample extractor MiniSampler, the sample liquid only comes into contact with the suction hose and the original sample bottle. Cross-contamination is largely avoided. The bottles are sealed immediately after sampling. The disposable suction hose can be replaced quickly. This ensures that unadulterated samples are taken each time. If necessary, a new tube is used for each sample. Due to its small diameter (8 mm), the tubing can also be guided to narrow, angled sample locations. Two versions: PE MiniSampler for general sampling, max. suction height 5 m\* PTFE MiniSampler for ultra-pure sampling, max. suction height 5 m\* \* depending on the medium Also suitable for viscous oils Complete in practical carrying case Not suitable for aggressive acids. Vapors could destroy pump.

# Additional Information

No.	KL-0641
Manufacturer (Brand)	Bürkle
VGKL number	100450110

