



Bürkle® SteriPlast Kit?sampling set, spatula+bag,sterile

Product Images



Description

The SteriPlast® Kit sampling set consists of ten sample spatulas or sample scoops and ten SteriBag Premium tamper-evident sample bags. The SteriPlast® Kit is thus the ideal combination for sampling where purity and sterility are particularly important. All parts are gamma-sterilized and are suitable for use with food. The scoops or spatulas are cleanroom-produced and individually packaged. The SteriBag Premium sample bag has an integrated seal and is liquid-tight. A unique tracking code both on the bag and on a tear-off field guarantees complete documentation and traceability. SteriPlast® Kit including scoop 50 ml or spatula 150 mm insertion depth, and SteriBag Premium 300 ml. Scoop 50 ml Bag SteriBag Premium 300 ml Spatula 150 mm insertion depth Bag SteriBag Premium 300 ml

Additional Information

No.	KL-0971
Manufacturer (Brand)	Bürkle
VGKL number	100458809

