



Marienfeld® HistoBond + blue

Product Images



Description

Slides HistoBond + blue, 100/pack HistoBond+ slides with adhesive surfaces for reliable attachment of tissue sections, glass, approx. 76x26x1 mm, with 90° ground edges, with blue writing field on one end, for applications according to EU Regulation 2017/746 for in vitro diagnostic medical devices (IVDR), with CE-mark, recommended use date and lot number for comprehensive information and traceability, packed per 100 pieces in plastic boxes, 2.000 pieces in outer carton

Additional Information

No.	MD-0211
Manufacturer (Brand)	Marienfeld

