



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Marienfeld® Haemacytometer Neubauer- improved light-lined

€87.70 As low

as

€74.70

**plus VAT &
Shipping**

Product Images



Description

Hemocytometer Neubauer-improved, bright lined made of optical special glass, with double net graduation and 0.1 mm depth, complete set in plastic case with 2 cover glasses, 1 blood mixing pipette each Thoma red and white, 2 color-coded silicone tubes of approx. 16 cm length, not registered for in-vitro diagnostic applications within the EU Please observe the guidelines of the accident prevention regulations when pipetting. We decline all responsibility when using the tubing for pipetting by mouth.

Additional Information

No.	MD-0121
Manufacturer (Brand)	Marienfeld
Transport temperature	Room temperature

